




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-130 USB Digital Microscope

Easy to use portable digital microscope with USB output interface
10X to 400X Magnification;

1.3 Mega Pixels CMOS sensor, interpolate to 2MP with software;

Manual Focus from 10mm to 500mm;

LED illumination with 8pcs LED lamps, brightness adjustable;

Photo management software, Auto White Balance and Exposure;

Powered by USB 2.0 cable, 5V DC;

Metal stand to hold the microscope steady;

Multi-language support: English, Germany, Spanish, French, Russian;

Operating System: Windows Win7, Win 8, Win 10 32bit&64 bit, Mac OS X 10.5 or above.