



LABOMED INC.[®]

Labomed Inc.

Mailing Address:

2728 S La Cienega Blvd. Los Angeles, CA 90034
U.S.A.

Emails: spectro@labomed.com,
labomedinc@gmail.com

Tel No.: (310) 202-0811

Fax No.: (310) 202-7286

[LB-122] Portable Monocular Inverted Biological Microscope



Portable Inverted Biological Microscope is designed to observe liquid specimens and other filed samples. It is portable, easy to operate and have very good performance.

[Download Product Brochure](https://labomedproducts.com/wp-content/uploads/2018/01/LB-122_USB_DIGITAL_MICROSCOPE_7-17-17_COMBINED-2.pdf) (https://labomedproducts.com/wp-content/uploads/2018/01/LB-122_USB_DIGITAL_MICROSCOPE_7-17-17_COMBINED-2.pdf)

[Get A Quote](https://www.labomedproducts.com/quote/?sku=[woo_sku]) ([https://www.labomedproducts.com/quote/?sku=\[woo_sku\]](https://www.labomedproducts.com/quote/?sku=[woo_sku]))

[\[print_link\]](#)

Features

Small portable microscope, simple to operate, good performance, high-definition.

Inverted design is suitable for observing liquid samples.

Powered by battery, with compact design, it is easy to carry out to inspect micro organism in water and factory effluent, outdoor observation of micro-plankton, pollen and spores.

Technical Specifications

Applications

LB-122 Portable Monocular Inverted Biological Microscope is mainly used to inspect micro organism in water and factory effluent, outdoor observation of micro-plankton, pollen and spores. So it can be used in water treatment factories, environmental protection department and also in biological laboratories.

Eyepiece Eyepiece 10×/18

Achromatic Objective 4× NA0.08/WD16

10× NA0.25/WD3.6

40× NA0.60/WD0.6

20× NA0.42/WD2.85

Nosepiece Triple Nosepiece

Stage Plain Stage with Slide Clips

Condenser&Illumination Bright field, LED Illumination(AA Battery X2)

Phase Contrast Observation, LED Illumination(AA Battery X2)

Moving Range 2.6mm

Leather Cover

Article printed from Labomed, Inc.: <https://labomedproducts.com>

URL to article: <https://labomedproducts.com/product/portable-monocular-inverted-biological-microscope/>

Copyright © 2017 Labomed Products. All rights reserved.